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INSTRUCTIONS FOR USE, STORAGE AND WARRANTY

Dear Customer,
Thank you very much for placing your trust in an Orliman product. Please read the instructions carefully. Keep these instructions and the packaging for future reference. If you have any questions or concerns, please contact your doctor, orthopaedic specialist or our customer service department.

ORLIMAN S.L.U. guarantees all its products as long as the original configuration has not been manipulated or altered except for the intended use as described in these instructions.

If the products are used in combination with other products, replacement parts or systems, make sure they are compatible and made by Orliman®. It does not guarantee any products with altered characteristics due to improper use, defects or breakage of any kind. To use the warranty, go to the place where you bought the product as they must complete the warranty information table on this instructions sheet. If any serious incidents related to the product occur, notify Orliman S.L.U. and the corresponding competent authority in your country.

Orliman would like to thank you for choosing this product and hopes you a speedy recovery.

REGULATIONS

[MD] This article is defined as an accessory of a class I medical device. A Risk Analysis (UNE EN ISO 14971) has been carried out, minimising the existing risks. Tests have been in accordance with European Regulation UNE-EN ISO 22523 on Prostheses and Orthoses.

INDICATIONS

- Allows the orthotic device, bandage or plaster to be used as intended (even while bathing or showering).
- Helps prevent infections and delayed recover due to contact with water.
- Prevents the need to replace a plaster due to damage caused by contact with water.

FITTING INSTRUCTIONS

For best therapeutic results considering different pathologies and to extend the useful life of the product, it is essential to choose the correct size for each patient or user. Excessive compression may be intolerable; adjusting the compression to be firm yet comfortable is recommended.

If the product must be adapted, any such adjustments must be done by an orthopaedic specialist or healthcare professional legally certified to do so who must make sure the end user or person responsible for fitting the product properly understands how it works and should be used.

When fitting the product, you must adhere to the following instructions:

When using the product on a child, fitting, use and removal must be supervised by an adult.

- Open the hole to insert the body segment to be covered.
- Carefully insert the limb with the orthotic device and/or plaster and slide the product over it until the area to be protected is covered. Try not to damage the material with fingernails or other sharp objects.
- Completely cover the orthotic device and/or plaster.
- Smooth all silicone folds so no liquid may enter and it becomes hermetically sealed around the skin above the orthotic device and/or plaster.

⚠ PRECAUTIONS

Before each use, check that all product components are present as per the fitting process. Periodically check the conditions of the product. If you observe any defect or anomaly, immediately report it to the issuing establishment.

This product is made of inflammable material. Do not expose the products to situations that could set them on fire. In the event of a fire, quickly get them off your body and use the proper resources to extinguish the fire.

For discomfort such as chafing, irritation and swelling, remove the product and see a doctor or orthopaedic specialist. The product should only be used on healthy skin. It is not recommended for use over open scars with swelling, redness or hotspots.

RECOMMENDATIONS-WARNINGS

The use of these products is conditioned by the indications. Although the product is not defined as a single-use device, using it on a single patient only is recommended and only for the intended purposes as described in these instructions or by a healthcare professional.

When disposing of the product and its packaging, you must strictly adhere to the legal regulations in your community.

This item is water-resistant. However, its use is not recommended for activities requiring immersion such as swimming and diving.

To avoid blood circulation problems, limit the length of use and make sure the product is not too tight.

Danger of asphyxia. Keep out of the reach of children. This product is not a toy.

Do not expose the product to the sun or other heat or flame sources.

Dry the top of the ring before removing the product.

Do not use sharp objects when handling such as earrings or rings.

If using with a plaster, make sure the rough edges are not in contact with the product as these edges may damage it.

Be particularly careful with moist, slipper and wet surfaces: do not walk or find support on them.

Do not reuse if the inside of the protector is dirty or damaged.

The plastic material may cause a feeling of heat.

RECOMMENDATIONS FOR STORAGE AND WASHING

When not using the product, store it in the original packaging in a dry place at room temperature. To dry the product, use a dry towel to absorb as much moisture as possible and let it dry at room temperature. Do not hang it up or iron the product and do not expose it to direct heat sources such as stoves, dryers, direct sun exposure, etc. When using or cleaning the product, do not use abrasive or corrosive substances, alcohol, ointments or liquid solvents. If not dried off properly, the detergent residue may irritate the skin and cause the product to deteriorate.

ES GARANTÍA 6 meses a partir de la compra			
EN GUARANTEE 6 months since purchase date			
FR GARANTIE 6 mois a partir de la date d'achat			
PT GARANTIA 6 meses a partir da compra			
ES Referencia	ES n° de lote	ES Fecha de compra	
EN Reference	EN lot no.	EN Purchase date	
FR Référence	FR lot non.	FR Date d'achat	
PT Referência	PT n° de lote	PT Data de compra	
ES Firma y sello de la ortopédia			
EN Stamp and signature of local stockist			
FR Cachet et signature du magasin spécialisé			
PT Nome e carimbo da ortopedia			
ES Para la validez de la garantía es necesario cumplimentar estos datos.		FR Pour valider la garantie, votre magasin spécialisé indiquer la date d'achat.	
EN To validate the warranty your local stockist should indicate the purchase date.		PT Para que a garantia seja válida e necessário completar com os dados.	

[EN] The manufacturing date is included in the batch number which can be found on the packaging label as **[LOT]**, in the following way: the second and third digits represent the year of manufacture and the fourth and fifth digits represent the month.

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